

# REAP-2 User Guide

## Table of Contents

<b>INPUT DATA.....</b>	<b>2</b>
DATASET INPUT REQUIREMENTS: .....	2
TRUNCATION STRATEGY: .....	2
EXAMPLE DATASET:.....	2
<b>PARAMETER SELECTION .....</b>	<b>2</b>
CHOOSING MODEL FEATURE .....	2
MODEL COMPARISON.....	3
PLOT SPECIFICS .....	4
DOWNLOAD REPORT .....	4

## Input data

### Dataset input requirements:

- The input dataset can be in csv, tsv and txt format
- The input dataset contains three columns: Concentration, Effect and Agent
- Columns in the input dataset should follow the order of Concentration, Effect and Agent

### Truncation strategy:

It is recommended that users normalize the response variable to the range of (0,1) by themselves. Otherwise, REAP will automatically truncate the values exceeding the boundaries to (0,1) using a truncation algorithm.

### Example dataset:

The example dataset can be downloaded through the following link:

[https://github.com/vivid225/REAP/blob/main/REAP/31780660\\_F1B\\_exempladata.csv](https://github.com/vivid225/REAP/blob/main/REAP/31780660_F1B_exempladata.csv)

Here is a screenshot of the dataset:

Auranofin	Cell viability	Agents
0.075	0.758	Mino
0.150	0.576	Mino
0.300	0.357	Mino
0.600	0.128	Mino
1.250	0.002	Mino
2.500	0.002	Mino

## Parameter selection

### Choosing model feature

#### 1. Log transform dose

The dose-response curve in REAP-2 is based on the median-effect equation. For statistical modeling, the dose-response curve is formulated as:

$$\log \frac{E}{1-E} = \beta_0 + \beta_1 \log d$$

where  $\beta_0$  and  $\beta_1$  are the intercept and slope parameters that determine a sigmoid dose-response relationship with respect to the effect  $E$  at the dose level  $d$ .

By unchecking the *log transform dose* option, the dose-response curve will be re-formulated as:

$$\log \frac{E}{1-E} = \beta_0 + \beta_1 d$$

## 2. Add potency estimation

The input value of the potency estimation is within (0, 100). By specifying the interested drug potency, REAP-2 will provide point estimation and standard deviation of the potency along with triangle signs specified in the dose-response curve plot (Red box in Fig. 1).




Figure 1 Dose-response curve plot. Triangles in the red box show IC50 values for each dose-response curve.

## Model comparison

### 1. Potency estimations

**Model Comparison**

Potency estimations

☒ yes

Slopes

☐ yes

By checking the *Potency estimations* option, REAP-2 will provide hypothesis testing on whether the **potency estimations** are the same for all the agents in the Summary tab.

### 2. Slopes

**Model Comparison**

Potency estimations


☐ yes

Slopes

☒ yes

By checking the *Slopes* option, REAP-2 will provide hypothesis testing on whether the **slope estimations** are the same for all the agents in the Summary tab.

### 3. Models



The image shows a dark-themed interface titled "Model Comparison". It contains two settings, each with a checked checkbox and the word "yes" to its right:

Setting	Value
Potency estimations	<input checked="" type="checkbox"/> yes
Slopes	<input checked="" type="checkbox"/> yes

By checking both the *Potency estimations* option and the *Slopes* option, REAP-2 will provide hypothesis testing on whether the **fitted models** are the same for all the agents in the Summary tab.

#### Plot specifics

1. Show values: points, means and StdDev  
In the Curve tab, users can choose what to include in the dose-response curve plot. *Points* represent the original raw data. *Means* provide mean estimations of the input dataset under each dose level. *StdDev* provide standard deviations of the response under each dose level. The default setting is a curve plot with mean and standard deviation of each dose level as in Fig. 1.
2. Log10(dose)  
*Log10(dose)* determines whether to log transform the dose level (x-axis) in the curve plot.
3. X-axis and y-axis names  
The default x-axis and y-axis names are based on the input dataset. Users can specify the names to change the axis labels of the curve plot.

#### Download report

There are three sections in the downloaded plot:

1. Setting includes all the parameter selections in REAP-2 for all the results in the report.
2. Dose-response curve plot saves the results under the Curve tab.
3. The tables in the report are saved from the results under the Summary tab.